

A Professional PerioCream® Product

PerioTabs®

**Helps reduce gingivitis, periodontitis,
peri-mucositis and peri-implantitis**

NITRADINE®

Patented Anti-Biofilm Formula



Made in Switzerland

A Professional PerioCream® Product

PerioTabs®



A BREAKTHROUGH IN ORAL CARE

What is PerioTabs®?

- A unique and easy to use brushing solution for gums and teeth.
- Helps reduce gingivitis and periodontitis.
- Helps reduce the risk of developing peri-mucositis and peri-implantitis.
- A safe, non-antibiotic biofilm removal formulation based on the world renowned NitrAdine®.
- No side effects such as teeth colouration, corrosion or burning sensations.



Product description

PerioTabs® helps reduce gingivitis, periodontitis, peri-mucositis and peri-implantitis through elimination of oral biofilm accumulation on gums, teeth and dental implants. The procedure consists of brushing the teeth and gums with the provided tablets dissolved in lukewarm water (Figure 1).

The formulation of PerioTabs® is based on the NitrAdine® anti-biofilm formula and is intended for use in the oral cavity on gums and teeth.

PerioTabs® contains 10 small effervescent tablets (1 tablet per day for 10 days). One tablet is to be dissolved in 15ml of lukewarm water in the provided container. Whilst the tablet is dissolving, the toothbrush should be immersed in the solution and left for 15 minutes. This will allow the toothbrush bristles to absorb the PerioTabs® solution (Figure 2) before brushing commences. The patient should then brush their teeth and gums with the solution for 2 minutes (Figure 3). This process should be repeated for 10 consecutive days.

PerioTabs® is a Medical Device Class I product, and patent pending in Europe and USA (application numbers EP15798353.7 and US 15/526,247).




Figure 1




Figure 2




Figure 3

What is NitrAdine®?

NitrAdine® is a unique tablet formulation. When dissolved in lukewarm water, it will result in a highly active solution that removes pathogenic biofilm from teeth, gums and implant materials to clinically insignificant levels.

The mode of action of NitrAdine®'s efficacy is based on the combination of a high concentration of surfactant (SLS) and other ingredients. The surfactant denatures the biofilm/micro-organism membrane proteins, whilst the other ingredients result in a slow release of a non-toxic concentration of only 0.02% hypochlorite, killing all micro-organisms (bacteria, fungi and viruses) present within the biofilm. (NitrAdine® is not hypochlorite. Hypochlorite, also called bleach, is a 2-10% concentration of NaOCl, which is toxic and cannot be used in the mouth).

Several independent academic *in-vitro* studies have demonstrated that NitrAdine® induces a significant reduction of biofilm formed by *Candida albicans*, *Staphylococcus aureus*, *Streptococcus mutans* and *Pseudomonas aeruginosa* on oral medical devices. In addition, an anti-viral activity against *Herpes simplex* was observed (Ref 1-8). Recent clinical studies in patients suffering from chronic periodontitis demonstrated that 10 days brushing with PerioTabs® induces a significant reduction of periodontopathogenic bacteria including *Porphyromonas gingivalis*, *Prevotella intermedia* and *Aggregatibacter actinomycetemcomitans* (Ref 9). NitrAdine® also reduces microorganisms present on toothbrush heads.

Additional properties of NitrAdine®:

- Non-toxic to gums (non-irritating, non-sensitizing)
- Non-antibiotic
- Non-corrosive to metal implants
- Non-erosive to dentine or enamel
- Compatible with all toothbrush heads




Efficacy of PerioTabs®

Clinical trial results


Multicentre, open label, randomised, clinical study to evaluate and compare the efficacy and safety of PerioCream®, a Periodontal Paste Dressing and Gum Brushing Solution (PerioTabs®) as an adjunct treatment to scaling and root planing (SRP) in patients suffering from chronic periodontitis.

Summary of Gingival Index progression over the treatment period




**Gingival Index was assessed using (Loe and Silness; 1963) index system for the assessment of the gingival condition and record qualitative changes on the gingival. It scores the marginal and interproximal tissues separately on the basis of 0 to 3.

Summary of Plaque Index progression over the treatment period




*Plaque levels were assessed using the Turkey modification of the Quigley-Hein Plaque Index.

Halitosis sensation after 11 days



Summary of Periodontal Bleeding Index progression over the treatment period




***Periodontal Bleeding Index was measured using Muhleman and Son 1971 score chart, which described the sulcus bleeding.

Microbial gum infection over the treatment period

Control group: 66.7% were infected at the start of the treatment. Out of these 66.7%, 100% were still infected after 11 days.

Treatment group: 60% were infected at the start of the treatment. Out of these 60%, only 25% were still infected after 11 days.



Efficacy of PerioTabs®

Case studies results



Before SRP



After SRP / 10 days
brushing with PerioTabs®



Post SRP



Post SRP / 10 days
brushing with PerioTabs®



Before



After 10 days brushing
with PerioTabs®




Before



After 5 days brushing
with PerioTabs®



Before



After 5 days brushing
with PerioTabs®

Efficacy of PerioTabs®

Case studies results



Before



Before



Before



After 10 days brushing
with PerioTabs®



After 10 days brushing
with PerioTabs®



After 10 days brushing
with PerioTabs®

Scientific references on NitrAdine®

1. Glass RT., Bullard JW., Conrad RS., Blewett EL. (2004). Evaluation of the sanitization effectiveness of a denture-cleaning product on dentures contaminated with known microbial flora. An in vitro study. *Quintessence Int* ; 35: 194–199.
2. Coenye T., De Prijck K., De Wever B., and Nelis H.J. (2008). Use of the modified Robbins device to study the in vitro biofilm removal efficacy of NitrAdine®, a novel disinfecting formula for the maintenance of oral medical devices. *J. Appl. Microbiol.* 105, 733–740.
3. Silva-Lovato C.H., De Wever B., Adriaens E., Paranhos H., Watanabe E., Pisani M.X., Souza R.F. de, and Ito I.Y. (2010). Clinical and antimicrobial efficacy of NitrAdine®-based disinfecting cleaning tablets in complete denture wearers. *J. Appl. Oral Sci.* 18, 560–565.
4. Vento-Zahra E., De Wever B., Decelis S., Mallia K., and Camilleri S. (2011). Randomized, double-blind, placebo-controlled trial to test the efficacy of NitrAdine® tablets in maxillary removable orthodontic appliance patients. *Quintessence Int.* 42, 37–43.
5. Glass R.T., Conrad R.S., Köhler G.A., Warren A.J., and Bullard J.W. (2011). Microbiota found in protective athletic mouthguards. *Sports Health* 3, 244–248.
6. Decelis S., Camilleri S., Zahra E.V., Scerri E., and De Wever B. (2012). The effect of NitrAdine® on the Candida levels of maxillary removable appliances. *Quintessence Int.* 43, 239–245.
7. Fathi H., Martiny H., Jost-Brinkmann P-G., (2015). Efficacy of cleaning tablets for removable orthodontic appliances. *J. Orofac. Orthop. / Fortschritte der Kieferorthopädie* 76, 143–151.
8. Coimbra F.C.T., Salles M.M., De Oliveira V.C., Macedo A.P., Da Silva C.H.L., Pagnano V.O., and Paranhos H. de F.O. (2016). Antimicrobial efficacy of complete denture cleansers. *Am. J. Dent.* 29, 149–153.
9. Sakly A., Jutla Y., Jutla B., Haddad I., De Wever B., Satia M., and Bogaert J.P. (2017). Clinical and microbiological effects of PerioCream® a periodontal dressing combined with an antimicrobial gum brushing solution after scaling and root planing : A multi-center randomized clinical trial. *J Dent Oral Care* 2, 1–6.

A Professional PerioCream® Product

PerioTabs®



Heiligkreuz 16
FL-9490 Vaduz
Liechtenstein, Europe

Phone: +423 232 78 15
Fax: +423 232 51 88
Email: administration@bonyf.com



www.periotabs.com

 **Made in Switzerland**